

General Surgical Instruments Catalogue



 Proud to be a family-run British company

**We are manufacturing
and supplying instruments
all over world.**



Surgical Holdings
Instruments



Symmetry
Surgical



Hermann
Laparoscopic



Bissinger
Bipolar



Light
Leads



DIN Baskets
& Containers



Sterile Packed
K-Wires



Custom-made
Instruments



**Working in close partnership with our customers,
we manufacture innovative and bespoke instrumentation.**

We have over 50 years of knowledge and experience in our workshop. Our team of highly skilled technicians undertake a minimum apprenticeship of 5 years and are trained to manufacture instruments, rather than just repair them. This gives us 360 degree understanding of the challenge and the importance of correctly functioning medical devices.

We care immensely about the work we do, never forgetting that instruments we create and repair are responsible for saving lives.

We insist on controlling the quality of every item we manufacture, from the materials used to the rigorous procedures we apply to testing. Our range of supporting products are serviced and repaired in-house, ensuring we are always able to exceed our customers' expectations.

Customers can receive training, participate in factory visits or simply call their local representative for guidance with instrumentation or management.

Skill-sharing has been at the heart of our business for four generations and we frequently arrange workshops and attend events at which we can educate others.

Quality Assurance

We are an ISO 13485 accredited company.

All Wellcare Surgical products undergo rigorous testing and inspection to applicable harmonised standards before leaving the factory. We offer CE marking on all Class I instrumentation, Class IIa and IIb implants, and can provide certificates of conformity for all materials used. Complete traceability throughout all aspects of the production process provides our customers with a fully documented quality control system.

All instrumentation carries a lifetime guarantee from the date of purchase for any fault found as a result of poor materials. This warranty is void, however, if repairs are carried out by anyone other than Wellcare Surgical repair personnel, or if it is deemed to have been misused.

Please contact us or visit our website for current reprocessing instructions.



wellcaresurgicals@live.co.uk


+44 7405903710 wellcaresurgicals.co.uk

General Instruments

Artery Forceps	4
Spring Forceps	7
Tissue Forceps	9
Sponge Holding and Sterilising Forceps	11
Scissors	12
Needle Holders	17
Combined Scissor and Needle Holders	21
Bag and Towel Clips	22
Retractors	23
Dissectors and Probes	31
Miscellaneous	32
Hollow Ware	33

Adson Artery Forceps

Box joint, 180mm long



Serrated jaws

Code	Description
WCS/034	Straight
WCS/034a	Curved




1-2 teeth

Code	Description
WCS/034b	Straight
WCS/034c	Curved

Cushing Artery Forceps

Box joint, fine, half serrated jaws, 140mm long




Code	Description
WCS/001	Straight
WCS/002	Curved



Cairns Artery Forceps

Box joint, half serrated jaws
145mm long



Serrated jaws

Code	Description
WCS/415	Straight




1-2 teeth

Code	Description
WCS/020	Curved
WCS/021	Curved to side

Dunhill Artery Forceps

Box joint, fine serrations



Straight


Code	Description
WCS/065	135mm long
WCS/065a	190mm long

Curved

Code	Description
WCS/043	135mm long
WCS/468	180mm long

Crile Artery Forceps

Box joint, serrated jaws



Straight

Code	Description
WCS/014	140mm long
WCS/014a	160mm long




Curved

Code	Description
WCS/015	140mm long
WCS/015a	160mm long

Lawson Tait Artery Forcep

Box joint, straight jaw




Code	Description
WCS/SH5924	140mm long

Artery Forceps

Halstead Mosquito Forceps

Standard pattern, haemostatic, box joint



Serrated jaw

Code	Description
WCS/007	Straight 125mm long
WCS/008	Curved 125mm long
WCS/009	Fully curved 125mm long

1-2 teeth

Code	Description
WCS/010	Straight 125mm long
WCS/012	Straight 190mm long
WCS/011	Curved 125mm long
WCS/013	Curved 190mm long

Micro, 100mm long

Code	Description
WCS/003	Straight, serrated
WCS/004	Curved, serrated
WCS/005	Straight, 1-2 teeth
WCS/006	Curved, 1-2 teeth

Kelly Artery Forceps

Box joint, half serrated jaws




Straight

Code	Description
WCS/022	140mm long
WCS/023	160mm long

Curved

Code	Description
WCS/024	140mm long
WCS/025	160mm long

Fraser Kelly Artery Forceps Curved




Straight

Code	Description
WCS/468a	180mm

Heiss Artery Forceps

Box joint, 205mm long



Serrated jaws


Code	Description
WCS/016	Straight
WCS/017	Curved

1-2 teeth

Code	Description
WCS/016a	Straight
WCS/017b	Curved

Kocher Artery Forceps

Box joint, 1-2 teeth



Straight

Code	Description
WCS/052	130mm long
WCS/053	150mm long
WCS/054	180mm long
WCS/055	200mm long
WCS/056	230mm long


Curved

Code	Description
WCS/057	150mm long
WCS/058	180mm long
WCS/059	200mm long
WCS/060	230mm long

Artery Forceps

Moynihan Artery Forceps

Box joint, half-serrated jaws



Straight


Code	Description
WCS/018	145mm long
WCS/068	180mm long

Curved

Code	Description
WCS/019	145mm long
WCS/069	180mm long

Rochester-Pean Artery Forceps

Box joint



Straight


Code	Description
WCS/044	165mm long
WCS/045	190mm long
WCS/046	205mm long
WCS/047	240mm long

Curved

Code	Description
WCS/048	165mm long
WCS/049	180mm long
WCS/050	205mm long
WCS/051	240mm long

Roberts Artery Forceps


Screw joint, curved on flat, serrated



Code	Description
WCS/070	230mm long

Spencer Wells Artery Forceps

Box joint, serrated jaws



Straight

Code	Description
WCS/035	127mm long
WCS/036	152mm long
WCS/037	177mm long
WCS/038	203mm long


Curved

Code	Description
WCS/039	127mm long
WCS/040	152mm long
WCS/041	177mm long
WCS/042	203mm long

OTHER SIZES AVAILABLE

Dandy Artery Forceps


Box joint, curved in plane



Code	Description
WCS/076	140mm long

Carmalt Angiotribe Forceps

Box joint, atraumatic, light pattern



Straight


Code	Description
WCS/062	165mm long
WCS/063	200mm long
WCS/064	230mm long

Curved

Code	Description
WCS/061	165mm long

Spring Forceps


Adsons Dissecting Forceps



127mm long

Code	Description
WCS/257	1-2 teeth
WCS/258	Lightweight, micro
WCS/259	Serrated
WCS/260	Serrated, lightweight, micro


Block End Dissecting Forceps



Serrated

Code	Description
WCS/255	127mm long
WCS/256	152mm long
WCS/256a	203mm long

Bonney Dissecting Forceps




178mm long

Code	Description
WCS/267	Serrated, 1-2 teeth
WCS/268	1-2 teeth

Charnley-McIndoe Dissecting Forceps


150mm long



Code	Description
WCS/ND45a	2-3 teeth

DeBakey Atraumatic Dissecting Forceps

Straight, atraumatic jaws



1.5mm

Code	Description
WCS/238a	150mm long
WCS/240a	240mm long

2.5mm

Code	Description
WCS/238b	150mm long
WCS/240b	240mm long

2mm


Code	Description
WCS/238	150mm long
WCS/241	180mm long

2.8mm

Code	Description
WCS/238c	150mm long
WCS/240c	240mm long
WCS/241c	180mm long

Dissecting Forceps

Turnover end, serrated jaws



Code	Description
WCS/242	127mm long
WCS/243	152mm long

Code	Description
WCS/244	178mm long
WCS/245	203mm long

Dunhill Dissecting Forceps


Light, fine serrated jaws



Code	Description
WCS/457	145mm long

Fredrick Coole Dissecting Forceps

Extra fine serrated jaws




125mm long

Code	Description
WCS/066	Straight
WCS/067	Curved

Spring Forceps

Gillies Dissecting Forceps


1-2 teeth



Code	Description
WCS/052	152mm long

McIndoe Dissecting Forceps


152mm long



Code	Description
WCS/476	Serrated
WCS/477	1-2 teeth

Jackson Burrows Dissecting Forceps


140mm long



Code	Description
WCS/N037	Tongue and groove, 1-2 teeth
WCS/N037a	Tongue and groove only

McIndoe Modified Dissecting Forceps


152mm long, 1mm tip, modified for fine powerful dissection with strong action tip



Code	Description
WCS/5961	Serrated


Jean Dissecting Forceps

Serrated jaws, 140mm long



Code	Description
WCS/269	Fine
WCS/270	1-2 teeth
WCS/271	2-3 teeth, 178mm long

Treves Dissecting Forceps



1-2 teeth


Code	Description
WCS/251	127mm long
WCS/252	152mm long

2-3 teeth

Code	Description
WCS/253	127mm long
WCS/254	152mm long

Lanes Dissecting Forceps


1-2 teeth on serrated jaws, 140mm long



Code	Description
WCS/246	127mm long
WCS/247	152mm long
WCS/248	178mm long

Code	Description
WCS/249	2-3 teeth, 178mm long
WCS/250	2-3 teeth 203mm long

Waugh Dissecting Forceps



Fine, serrated

Code	Description
WCS/261	152mm long
WCS/262	178mm long
WCS/263	203mm long


Fine, 1-2 teeth

Code	Description
WCS/264	152mm long
WCS/265	178mm long
WCS/266	203mm long

Tissue Forceps

Allis Tissue Forceps

Box joint



150mm long


Code	Description
WCS/272	3-4 teeth
WCS/273	4-5 teeth

190mm long

Code	Description
WCS/414	5-6 teeth

Duval Tissue Forceps


Box joint



Code	Description
WCS/279	13mm wide, 203mm long

Littlewood Tissue Forceps


Box joint, 2-3 teeth



Code	Description
WCS/275	190mm long

Lane Tissue Forceps


Box joint, 1-2 teeth, fenestrated



Code	Description
WCS/280	127mm long
WCS/281	152mm long
WCS/282	197mm long

Babcock Tissue Forceps


Box joint, serrated jaws



Code	Description
WCS/276	160mm long
WCS/277	200mm long
WCS/278	230mm long

Rutherford Morrison Tissue Forceps

Box joint




Code	Description
WCS/283	4-5 teeth, 152mm long

Tissue Forceps

Shaw Tissue Forceps

Box joint




Code	Description
WCS/458	1-2 teeth, 178mm long



Guys Gallabin Suture-Holding Forceps


Screw joint, angled flat blades



Code	Description
WCS/450	125mm long
WCS/4402	150mm long

Russian Dissecting Forceps

With fenestrated end



Code	Description
WCS/SH-5503b	150mm long
WCS/SH-5503	200mm long
WCS/403	230mm long
WCS/SH-5503a	250mm long



Special Offer


on SSD Essentials Pack

Top up your trays with these essentials...

CONTACT US FOR OUR OFFER PRICE


Halstead Mosquito Forceps

Straight/Curved 125mm (5")




Ramsey Sponge Holding Forceps

Straight, Serrated Jaws
241mm
(9½ ")




Scalpel Blade

Handle
No. 3




Spring Clip

Serrated Pattern,
Stainless Steel for
Fixing Slings to Leg
Splints, also for Towel
Clipping
6.3cm (2½ ") Wide



Spencer Wells Artery Forceps

Straight, Serrated
Jaws




Email us to place your order
wellcaresurgicals@live.co.uk

Sponge Holding and Sterilising Forceps

Rampley Sponge Holding Forceps


Box joint, straight, serrated jaws




Code	Description
WCS/353	178mm long
WCS/354	241mm long

Lister Sinus Forceps

Box joint, serrated tip




Code	Description
WCS/367	127mm long
WCS/368	152mm long
WCS/369	178mm long




Foerster Sponge Holding Forceps

Box joint, straight, serrated jaws




Code	Description
WCS/354a	178mm long



Cheatle Sterilizer Forceps

Screw joint




Code	Description
WCS/364	267mm long

Harrison Bowl Steriliser Forceps


Screw joint, angled, double jaw

Code	Description
WCS/365	356mm long
WCS/366	457mm long



Premium Quality Stainless Steel Scissors

Dressing Scissors



Blunt/Blunt

Straight

Code	Description
WCS/133	130mm long
WCS/136	146mm long
WCS/139	165mm long

Curved

Code	Description
WCS/142	130mm long
WCS/143	146mm long
WCS/144	165mm long
WCS/145	185mm long
WCS/404	147mm long

With one serrated blade

Blunt/Sharp

Straight

Code	Description
WCS/134	130mm long
WCS/137	146mm long
WCS/140	165mm long

Sharp/Sharp


Straight

Code	Description
WCS/135	130mm long
WCS/138	146mm long
WCS/141	165mm long

Curved

Code	Description
WCS/146	130mm long

Metzenbaum Scissors



Straight

Code	Description
WCS/405	180mm long
WCS/406	147mm long
WCS/406a	200mm long
WCS/406b	230mm long
WCS/406c	250mm long
WCS/406d	280mm long

Curved

Code	Description
WCS/407	185mm long
WCS/408	147mm long
WCS/408a	200mm long
WCS/408b	230mm long
WCS/408c	250mm long
WCS/408d	280mm long

Mayo Scissors



Straight

Code	Description
WCS/409	140mm long
WCS/409a	160mm long
WCS/411	190mm long
WCS/411a	230mm long

Curved


Code	Description
WCS/410	140mm long
WCS/410a	160mm long
WCS/412	190mm long
WCS/412a	230mm long

Dressing Scissors, Left-Handed

With one black bow, 145mm long

Blunt/Blunt

Code	Description
WCS/192	Straight
WCS/193	Curved



Blunt/Sharp

Code	Description
WCS/194	Straight
WCS/195	Curved


Sharp/Sharp

Code	Description
WCS/196	Straight
WCS/197	Curved

Littler Dissecting Scissors

Fine pattern, with eye for suture

Code	Description
WCS/001f	115mm long



Kilner Dissecting Scissors


Fine pattern

Code	Description
WCS/001a	115mm long



Fergusson Scissors


Angled on flat



Code	Description
WCS/148	180mm long

Paynes Suture Scissors


Code	Description
WCS/152	100mm long



Ligature Stitch Scissors


Fine pattern

Code	Description
WCS/153	115mm long



Scissors


Lister Bandage Scissors



Code	Description
WCS/154	147mm long
WCS/155	178mm long

Lister Bandage Scissors


Angled, with probe end



Code	Description
WCS/361	147mm long
WCS/362	178mm long

McIndoe Scissors

178mm long




Code	Description
WCS/149	Straight
WCS/150	Curved on flt at
WCS/150sc	Supercut, curved




Thomas Strabismus Scissors

114mm long




Code	Description
WCS/030	Straight
WCS/029	Curved

Heath Stitch Scissors




Code	Description
WCS/151	152mm long

“ Careless” Scissors



Code	Description
WCS/474	165mm long

Tubing Scissors



Code	Description
WCS/147	130mm long

Utility Scissors




Code	Description
WCS/466	190mm long

Tungsten Carbide Hard Insert Scissors

Dressing Scissors

Tungsten carbide hard insert,
blunt/blunt



Straight


Code	Description
WCS/430	130mm long
WCS/432	145mm long
WCS/434	165mm long
WCS/436	185mm long

Curved

Code	Description
WCS/431	130mm long
WCS/433	145mm long
WCS/435	165mm long

Mayo Dressing Scissors

Tungsten carbide hard insert



Straight


Code	Description
WCS/419	140mm long
WCS/419a	160mm long
WCS/421	190mm long
WCS/421a	230mm long

Curved

Code	Description
WCS/420	140mm long
WCS/420a	160mm long
WCS/422	190mm long
WCS/422a	230mm long

Dressing Scissors

Tungsten carbide hard insert,
blunt/sharp



Straight


Code	Description
WCS/423	130mm long
WCS/425	145mm long
WCS/427	165mm long
WCS/429	185mm long

Curved

Code	Description
WCS/424	130mm long
WCS/426	145mm long
WCS/428	165mm long

Mayo-Harrington Dressing Scissors

Tungsten carbide hard insert



230mm long


Code	Description
WCS/417	Straight
WCS/418	Curved

200mm long

Code	Description
WCS/162a	Straight
WCS/413b	Curved

Dressing Scissors

Tungsten carbide hard insert sharp/sharp



Straight


Code	Description
WCS/437	130mm long
WCS/439	145mm long
WCS/441	165mm long
WCS/443	185mm long

Curved

Code	Description
WCS/438	130mm long
WCS/440	145mm long
WCS/442	165mm long

Metzenbaum Scissors

Tungsten carbide hard insert



230mm long


Code	Description
WCS/445	145mm long
WCS/448	185mm long
WCS/471a	230mm long
WCS/471b	250mm long
WCS/471c	280mm long

200mm long

Code	Description
WCS/446	145mm long
WCS/449	185mm long
WCS/449a	200mm long
WCS/449b	230mm long
WCS/449c	250mm long
WCS/449d	280mm long

Baby Metzenbaum Scissors


Tungsten carbide,
hard insert, 145mm long



Code	Description
WCS/445a	Straight
WCS/446a	Curved

Ligature Scissors

Tungsten carbide
hard insert, curved




Code	Description
WCS/444	145mm long
WCS/447	185mm long

Tungsten Carbide Supercut Scissors

Mayo Harrington Scissors

Supercut, tungsten carbide



230mm long


Code	Description
WCS/162	Straight
WCS/413	Curved

200mm long

Code	Description
WCS/162b	Straight
WCS/413b	Curved

Mayo Scissors

Supercut, tungsten carbide



Straight


Code	Description
WCS/163	140mm long
WCS/163a	160mm long
WCS/165	190mm long
WCS/165a	230mm long

Curved

Code	Description
WCS/164	140mm long
WCS/164a	160mm long
WCS/166	190mm long
WCS/166a	230mm long

Metzenbaum Dissecting Scissors

Supercut, tungsten carbide



Straight


Code	Description
WCS/188tc	185mm long
WCS/188tca	200mm long
WCS/188tcb	230mm long
WCS/188tcc	250mm long
WCS/188tcd	280mm long

Curved

Code	Description
WCS/191tc	147mm
WCS/189tc	185mm
WCS/186tca	200mm
WCS/188tcb	230mm
WCS/188tcc	250mm
WCS/188tcd	280mm

Dressing Scissors

Supercut, tungsten carbide



Blunt/Blunt

Straight

Code	Description
WCS/174	130mm long
WCS/176	145mm long
WCS/178	165mm long
WCS/180	185mm long

Curved

Code	Description
WCS/175	130mm long
WCS/177	145mm long
WCS/179	165mm long
WCS/180a	185mm long

Blunt/Sharp

Straight

Code	Description
WCS/167	130mm long
WCS/169	145mm long
WCS/171	165mm long
WCS/173	185mm long

Curved

Code	Description
WCS/168	130mm long
WCS/170	145mm long
WCS/172	165mm long
WCS/173a	185mm long

Sharp/Sharp

Straight

Code	Description
WCS/181	130mm long
WCS/183	145mm long
WCS/185	165mm long
WCS/187	185mm long

Curved

Code	Description
WCS/182	130mm long
WCS/184	145mm long
WCS/186	165mm long
WCS/187a	185mm long

Guide for Gold and Black Bows

It's important to know the difference for the intended use of each scissor:



2 Gold Bows
Tungsten carbide insert or padding for a scissor



2 Black Bows
Supercut edge



1 Black and 1 Gold Bow
Supercut tungsten carbide



1 Black Bow (left ring)
Left-handed instrument

Supercut Operating Scissors



WITH TWO BLACK BOWS
FOR IDENTIFICATION

Dressing Scissors

Supercut, blunt/blunt




145mm long

Code	Description
WCS/160	Straight
WCS/161	Curved

Dressing Scissors

Supercut, blunt/sharp



145mm long

Code	Description
WCS/156	Straight
WCS/157	Curved

Dressing Scissors

Supercut, sharp/sharp




145mm long

Code	Description
WCS/158	Straight
WCS/159	Curved

Mayo Scissors

Supercut



Straight


Code	Description
WCS/166a	140mm long
WCS/166e	160mm long
WCS/166b	170mm long

Curved

Code	Description
WCS/166c	140mm long
WCS/166f	160mm long
WCS/166d	170mm long

Boyd Scissors


Supercut, blunt/blunt



Code	Description
WCS/047	180mm long

Metzenbaum Dissecting Scissors

Supercut



Straight (length)

Code	Description
WCS/190	147mm
WCS/188	185mm
WCS/188a	200mm
WCS/188b	230mm
WCS/188c	250mm
WCS/188d	280mm

Curved (length)

Code	Description
WCS/191	147mm
WCS/189	185mm
WCS/186a	200mm
WCS/186b	230mm
WCS/186c	250mm
WCS/186d	280mm

Stevens Tenotomy Scissors

Supercut



Curved/blunt

Code	Description
WCS/424	114mm long
WCS/424a	125mm long
WCS/424gb	135mm long

Curved/sharp

Code	Description
WCS/425	114mm long
WCS/425a	140mm long
WCS/425b	160mm long

With cutout, 110mm long

Code	Description
WCS/426	Straight/blunt, 114mm long
WCS/427	Straight/sharp, 114mm long
WCS/427c	Straight/sharp, 180mm long

Code	Description
WCS/427a	Straight/blunt
WCS/427b	Curved/blunt

Strabismus Scissors

Supercut




114mm long

Code	Description
WCS/475	Straight
WCS/475a	Curved

Iris Scissors

Supercut




114mm long

Code	Description
WCS/476	Straight
WCS/476a	Curved

Needle Holders

Bruce Clarke Needle Holder


Screw joint



Code	Description
WCS/205	127mm long

Lawrence Needle Holder


Fine jaw, 150mm long



Code	Description
WCS/012a	
WCS/012b	Tungsten carbide jaw

Crile-Wood Needle Holder


Straight, baby, fine



Code	Description
WCS/199	150mm long

Needle Holder


Surgical Holdings pattern, box joint, lightweight, with curved plain jaws and one o set bow for thumb



Code	Description
WCS/209	152mm long

Kilner Needle Holder


Box joint, cross-serrated jaw



Code	Description
WCS/200	165mm long
WCS/201	133mm long

MacPhail Needle Holder

Screw joint, with one copper jaw




Code	Description
WCS/206	159mm long
WCS/207	178mm long
WCS/208	209mm long

Needle Holders

Mathieu Needle Holder


Screw joint



Code	Description
WCS/210	172mm long
WCS/211	209mm long
WCS/212	243mm long


Naughton Morgan Needle Holder

Box joints, cross serrated jaws



Mayo Needle Holder

Box joint, straight




Tungsten carbide

Code	Description
WCS/399	185mm long
WCS/400	254mm long

Code	Description
WCS/401	185mm long
WCS/402	254mm long

Rochester-Mayo Needle Holder

Screw joint



Wide cross serrated jaws

Code	Description
WCS/224	159mm long
WCS/398	203mm long

Narrow cross serrated jaws

Code	Description
WCS/225	159mm long


Code	Description
WCS/213	178mm long

Choosing the most suitable needle size

Our tungsten carbide inserts are available in a variety of styles. Choosing the most suitable needle size for each procedure helps give the surgeon maximum grip and control.


Coarse 0.65mm

For powerful needle holder patterns and wire twisters. Maximum holding capabilities for heavy needles or wire.




Micro 0.2-0.3mm

For micro needle holders such as Castroviejo, Microvascular and Mini-Ryder. For delicate needles and sutures between 7-0 and 10-0.



Standard 0.5mm

Used, for example, in Mayo-Hegar, Mathieu Needle Holders. For needles 3-0 and larger.




Smooth

Smooth tungsten carbide surface reduces damage to delicate needles and suture of 6-0 and finer.



Fine 0.4mm

For delicate needle holders such as Crile-Wood and DeBakey. For needles and sutures between 4-0 and 6-0.




Tungsten Carbide Needle Holders

 WITH TWO GOLD BOWS FOR IDENTIFICATION

Berry Needle Holder/Wire Twister


Tungsten carbide jaws



Code	Description
WCS/467	200mm long

Finochietto Needle Holder


Tungsten carbide jaws



Code	Description
WCS/222	215mm long
WCS/223	265mm long

Crile-Wood Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



Code	Description
WCS/217	159mm long
WCS/218	178mm long
WCS/198	150mm long

Halsey Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



Code	Description
WCS/220	127mm long

Derf Needle Holder


Tungsten carbide insert



Code	Description
WCS/BUPHF002	115mm long

Kilner Needle Holder

Tungsten carbide insert




Code	Description
WCS/221	140mm long

Tungsten Carbide Needle Holders

 WITH TWO GOLD BOWS FOR IDENTIFICATION

Masing Needle Holder


Tungsten carbide insert, one angled bow



Code	Description
WCS/009R	130mm long

Nievert Needle Holder


Tungsten carbide insert, box joint, straight



Code	Description
WCS/219	127mm long

Ryder Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



0.2mm jaw pitch


Code	Description
WCS/SH-5197a	135mm long
WCS/SH-5197b	155mm long
WCS/SH-5197c	175mm long
WCS/SH-5197d	210mm long

0.4mm jaw pitch

Code	Description
WCS/SH-5198a	135mm long
WCS/SH-5198b	155mm long
WCS/SH-5198c	175mm long
WCS/SH-5198d	210mm long

DeBaKey Needle Holder


250mm long



Code
WCS/SC-5642a

Sarot Needle Holder


Diamond jaw, tungsten carbide, gold plated bows



Code	Description
WCS/451	181mm long

Tungsten Carbide Needle Holders


Features and Benefits



- Jaw**
Tungsten carbide inserts bonded into the jaws for a long lasting, secure grip
- Bevelled Edges**
Prevents suture from snagging
- Laser Etching**
Date and catalogue number are etched for easy re-ordering. Laser etch will not harm instrument like traditional methods
- Satin Finish**
Reduces glare
- Ratchets**
Precisely cut for effortless locking/unlocking
- Teeth**
Fine, hard teeth, diamond cut into an insert of hardened tungsten carbide for superior needle grasping and suture tying capabilities
- Vacuum Brazing**
Combines bonding of the tungsten carbide insert with the simultaneous vacuum hardening of the instrument for greater durability and longevity
- Shanks**
Precision balanced for comfort and control
- Gold Handles**
For easy identification

Mayo-Hegar Needle Holder

Tungsten carbide insert, box joint, satin finish, straight




Code	Description
WCS/472	130mm long
WCS/214	159mm long
WCS/215	178mm long
WCS/216	209mm long

LONGER SIZES AVAILABLE

Combined Scissor and Needle Holders

Olsen-Hegar Scissor/Needle Holder


With tungsten carbide insert



Code	Description
WCS/229	147mm long
WCS/230	165mm long

Olsen-Hegar Scissor/Needle Holder


Screw joint



Code	Description
WCS/226	140mm long
WCS/227	165mm long
WCS/228	190mm long

Gillies Scissor/Needle Holder

Right-hand, with tungsten carbide insert



Code	Description
WCS/235	158mm long

Gillies Scissor/Needle Holder


Heavy pattern, 158mm long



Code	Description
WCS/233	Right-hand
WCS/234	Left-hand

Gillies Scissor/Needle Holder


158mm long



Code	Description
WCS/231	Right-hand
WCS/232	Left-hand

Foster Gillies Scissor/Needle Holder

Fine pattern




Code	Description
WCS/236	127mm long

Bag and Towel Clips

Spring Clip

Serrated pattern, stainless steel, for fixing slings to leg splints, also for towel clipping



Code	Description
WCS/344	63mm wide

Spring Bag Clip


Serrated



Code	Description
WCS/345	26mm wide

Jones Towel Clip


Superior quality spring forcep type



Code	Description
WCS/346	89mm long

Shardles Towel Clip


Spring curved, cross action



Code	Description
WCS/347	77mm long

Backhaus Towel Clip


Superior quality, box joint



Code	Description
WCS/348	89mm long
WCS/348a	102mm long
WCS/349	130mm long

Moynihan Towel Clip


Screw joint, 2-2 teeth, for sides of incision



Code	Description
WCS/470	191mm long

Ball and Socket Towel Forceps

Box joint




Code	Description
WCS/350	89mm long
WCS/351	115mm long
WCS/352	140mm long

Retractors

Goligher Retractor

Complete with 2 centre blades:
 1 x 90mm wide x 50mm deep
 1 x 90mm wide x 90mm deep and 2
 pairs of lateral blades:
 1 pair 25mm wide x 45mm deep 1 pair
 50mm wide x 90mm deep



Code
WCS/305

Balfour Abdominal Retractor

Trivalve, frame 203mm

Child


102mm opening, 51mm deep lateral blades, 38mm
 deep centre blade

Adult

178mm opening, 64mm deep lateral blades, 44mm
 deep centre blade

Extra large


178mm opening, 102mm deep lateral blades, 108mm
 deep centre blade



Code	Description
WCS/373	Child
WCS/374	Adult
WCS/375	Extra large

Miles Retractor

Complete with 2 centre blades:
 1 x 90mm wide x 50mm deep
 1 x 90mm wide x 75mm deep and 2
 pairs of lateral blades:
 1 pair 80mm wide x 35mm deep 1 pair
 80mm wide x 75mm deep




Code	Description
WCS/376	340mm long

OTHER BLADE SIZES AVAILABLE

Balfour Abdominal Retractor


Bivalve, baby



Code	Description
WCS/372	Fully open 89mm

Berkeley Abdominal Retractor

Parallel action




Code	Description
WCS/461	Blades 35mm x 60mm deep, 215mm long

Turner Warwick Retractors

Turner Warwick Universal Self Retaining Retractor

Code	Description
WCS/520	With frames A/B/C/D and blades A/B/C/D/E/F/K/L
WCS/521-A	Circular frame A only 178mm
WCS/521-B	Circular frame B only 229mm
WCS/521-C	Circular frame C only 305mm
WCS/521-D	Oval frame D only 475mm x 330mm
WCS/521-E	Oval frame E only 370mm x 300mm



Turner Warwick Needleholder with Tungsten Carbide Jaws

Curved shanks with one short bow


Code	Description
WCS-5835	230mm long
WCS-5835A	200mm long

Blade A




Code	Description
WCS/522-A	Blade only - 19mm x 76mm

Blade B




Code	Description
WCS/522-B	Blade only - 38mm x 102mm

Blade C



Code	Description
WCS/522-C	Blade only - 57mm x 89mm


Blade D



Code	Description
WCS/522-D	Blade only - 57mm x 114mm


Turner Warwick Retractors

Blade E




Code	Description
WCS/522-E	Blade only - curved stainless steel blade 51mm x 38mm

Blade F




Code	Description
WCS/522-F	Blade only - curved stainless steel blade 76mm x 60mm

Blade G




Code	Description
WCS/522-G	Blade only - flat stainless steel blade 89mm x 60mm

Blade J




Code	Description
WCS/522-J	Blade only - single hook, 13mm deep

Blade K



Code	Description
WCS/522-K	Blade only - flat with clip 55mm x 90mm

Blade L




Code	Description
WCS/522-L	Blade only - flat with clip 55mm x 120mm

Retractors

Collin Retractor


110mm fully open, 235mm long



Code	Description
WCS/322	Complete with centre blade 55 x 40mm and one pair of blades 38mm x 60mm
WCS/322a	Without centre blade

Gelpi Retractor


Self retaining



Code	Description
WCS/O13a	140mm long
WCS/O13b	178mm long


Joll Thyroid Retractor

Self retaining, with lock-in nut




Code	Description
WCS/363	152mm long

Gossett Abdominal Retractor




Code	Description
WCS/303	Large, fully open 185mm long
WCS/304	Small, fully open 95mm long

Pozzi Abdominal Retractors



Code	Description
WCS/370	Swivel blades: 38 x 60mm, bivalve, 254mm long




Code	Description
WCS/371	Swivel blades: 38 x 59mm, trivalve, 203mm long

Retractors

Bowlby Dilator


Curved



Code	Description
WCS/495	145mm long

Travers Retractor


Self retaining with internal rack, 210mm long



Code	Description
WCS/323	4 x 5 teeth, 95mm fully open

Irwin Retractor


Self retaining, with cam rack, 4 x 4 teeth



Code	Description
WCS/324	178mm long

West Retractor


Self retaining, 140mm long



Code	Description
WCS/326	3 x 4 teeth, 60mm fully open
WCS/326cr	3 x 4 teeth, 60mm fully open, with cam rack

Norfolk and Norwich Retractor


4 x 5 teeth 50mm blades x 223mm long



Code	Description
WCS/325	Standard rack

Weitlander Retractor

3-4 teeth



135mm long


Code	Description
WCS/496	Blunt
WCS/496a	Sharp
WCS/496b	Sharp, 2-3 teeth
WCS/496c	Blunt, 2-3 teeth

165mm long

Code	Description
WCS/497	Blunt
WCS/497a	Sharp

Retractors

Kocher Retractor



230mm long


Code	Description
WCS/378	36 x 20mm blade

250mm long

Code	Description
WCS/378a	67 x 25mm blade
WCS/378b	80 x 42mm blade
WCS/378c	80 x 50mm blade
WCS/378d	80 x 60mm blade

Canny Ryall Navy Retractor


Double-ended, with one single hook blade 5mm wide, and one fenestrated blade



Code	Description
WCS/343	Fenestrated blade 16mm wide, 190mm long
WCS/343a	Fenestrated blade 19mm wide, 190mm long
WCS/343b	Fenestrated blade 25mm wide, 220mm long
WCS/343c	Fenestrated blade 30mm wide, 220mm long

Kocher Langenbeck Retractor

215mm long



Code	Description
WCS/336sp	25 x 6mm blade
WCS/379a	35 x 8mm blade
WCS/379b	35 x 11mm blade
WCS/379c	35 x 15mm blade
WCS/357sp	41 x 11mm blade

Code	Description
WCS/cutsp015	55 x 11mm blade
WCS/379d	70 x 15mm blade
WCS/379e	80 x 12mm blade
WCS/379f	80 x 16mm blade

Canny Ryall Army Retractor

Double prong, 13mm wide, 165mm long



Code	Description
WCS/342	Large, 19mm wide
WCS/342a	Small, 11mm wide



Langenbeck Retractor


Single ended



Code	Description
WCS/336	Blade 22mm x 6mm
WCS/337	Blade 45mm x 13mm
WCS/337a	Blade 65mm x 13mm
WCS/338a	Blade 75mm x 25mm
WCS/338a1	Blade 65mm x 25mm


Czerny Retractor

Double ended



Code	Description
WCS/341	171mm long

Landon Retractor




Code	Description
WCS/494	230mm



Lane Retractor

Double ended blade, 229mm long




Code	Description
WCS/339	13mm x 19mm wide
WCS/340	25mm x 32mm wide

Retractors

Morris Retractor


Single blade, 229mm long



Code	Description
WCS/330	38 x 51mm
WCS/331	51 x 51mm
WCS/332	64 x 51mm


Morris Retractor

Double ended, 216mm long




Code	Description
WCS/327	32 x 51mm and 38 x 51mm
WCS/328	44 x 51mm and 51 x 51mm
WCS/329	64 x 51mm and 38 x 51mm

Moynihan Retractor




Code	Description
WCS/462	Single ended, blade 13 x 22mm, 215mm long

Deaver Abdominal Retractors



Code	Description
WCS/306	Set of 5, 315mm overall length
WCS/307	25mm wide
WCS/308	25mm wide
WCS/309	38mm wide
WCS/310	51mm wide
WCS/311	76mm wide

Elliot (Modified Deaver) Retractor



Code	Description
WCS/308M	With fl at handle

Dyball and Lock Retractor

For liver




Code	Description
WCS/335	64 x 144mm, 241mm long

OTHER SIZES AVAILABLE - PRICE ON APPLICATION

Deaver Abdominal Retractors


With round handles



Code	Description
WCS/312	Set of 5, 305mm overall length
WCS/313	25mm wide
WCS/314	25mm wide
WCS/315	38mm wide
WCS/316	51mm wide
WCS/317	76mm wide

Kelly Retractor Blade

267mm long




Code	Description
WCS/333	38 x 191mm
WCS/334	57 x 191mm

Retractors

Durham Retractor


With fenestrated blade, 216mm long



Code	Description
WCS/300	10mm wide
WCS/301	14mm wide
WCS/302	20mm wide

Ballance (Volkman) Retractor

203mm long



Sharp

Code	Description
WCS/463	1 prong
WCS/463a	2 prong
WCS/463b	3 prong
WCS/463c	4 prong
WCS/463d	6 prong

Blunt

Code	Description
WCS/463e	1 prong
WCS/463f	2 prong
WCS/463g	3 prong
WCS/463h	4 prong
WCS/463i	6 prong

Volkman Retractor



Code	Description
WCS/400	Sharp, 220mm long
WCS/400a	Blunt, 220mm long
WCS/400b	Sharp, 245mm long

Kilner Retractor

Double ended, cat's paw



Code	Description
WCS/056	152mm long

Kocher Retractor



Code	Description
WCS/N036c	Sharp, 205mm long
WCS/N036	Sharp, 220mm long
WCS/N036b	Blunt, 220mm long

Senn-Miller Retractor

170mm long



Code
WCS/5881

Single Hook Retractor

Sharp



Code	Description
WCS/SH5671	152mm long
WCS/065b	Half size, fine

Gillies Skin Hook

Solid handle, 190mm long



Code	Description
WCS/065a	Standard, fine
WCS/065b	Half size, fine

Dissectors and Probes

MacDonald Dissector Double Ended



Code	Description
WCS/285	191mm long

Probe, Sinus Silver



Code	Description
WCS/464	152mm long

Brodie Director Probe Ended, Silver Plated



Code	Description
WCS/286	165mm long

Probe, Fish Tail Silver Plated



Code	Description
WCS/293	127mm long
WCS/294	152mm long
WCS/295	203mm long
WCS/296	254mm long

Watson Cheyne Dissector

Fine, double-ended, with probe end



Code	Description
WCS/287	127mm long
WCS/288	191mm long
WCS/288a	230mm long
WCS/288b	255mm long

Jobson Horne Probe Dressing (Aural)

With one flat tapered end, serrated edge and one small ring end, stainless steel



Code	Description
WCS/297	140mm long
WCS/298	178mm long

Probe, Malleable Silver Plated



Code	Description
WCS/289a	127mm long
WCS/289	152mm long
WCS/290	178mm long
WCS/291	203mm long
WCS/292	255mm long
WCS/292a	305mm long

Bookwalter Retractor Kit III

A table-fixed retractor system designed by surgeons to efficiently and reliably hold surgical instruments in place during surgery while allowing greater visibility.

Kit includes Bookwalter® Kit:
Table Post and Coupling, Horizontal Flex Bar, Segmented Rings, Tilt Ratchet Mechanism, Kelly Retractors, Blades and Ratchets

Contact us for full details at wellcaresurgicals@live.co.uk

Miscellaneous


Aneurysm Needle

Small, for infusions



Code	Description
WCS/472a	146mm long

Mayo Instrument Pin



Code	Description
WCS/355	102mm long
WCS/356	152mm long

Syme Aneurysm Needle



Code	Description
WCS/473	165mm long


Baron Scalpel Handle

130mm long



Code	Description
WCS/003	No 3
WCS/004	No 4

Ring Cutter



Code	Description
WCS/452	Circular saw blade, 165mm long
WCS/452a	Spare blade

Scalpel Handle

127mm long



Code	Description
WCS/007	No 3
WCS/008	No 4
WCS008aNo 4 L	
WCS/001aNo 7	

Endoscope Re-fibre

Wellcare surgical manufacture their own fibre bundles in house. We cover small and larger diameter scopes and only use high quality German fibre optic material and stainless steel tubing.

Quality components, globally sourced:


- ✓ Medical grade stainless steel tubing
- ✓ German engineered new glass fibre
- ✓ New fibre cone as standard
- ✓ New solderable sapphire distal window
- ✓ New light post colour ring
- ✓ New objective in 95% of cases
- ✓ Replacement 'O' ring to prevent leakage

Contact us with your requirements
wellcaresurgicals@live.co.uk

Hollow Ware

Kidney Dish


Stainless steel



Code	Description
WCS/380	152mm
WCS/381	203mm
WCS/382	254mm
WCS/383	305mm

Lotion Bowl


Without lid, stainless steel



Code	Description
WCS/384	102mm
WCS/385	152mm
WCS/386	203mm
WCS/387	254mm
WCS/388	305mm

Instrument Tray

Stainless steel



Code	Description
WCS/389	254 x 203 x 50mm
WCS/390	Lid only for 254 x 203 x 50mm
WCS/391	305 x 254 x 50mm
WCS/392	Lid only for 305 x 254 x 50mm
WCS/393	356 x 305 x 50mm
WCS/394	Lid only for 356 x 305 x 50mm

Cheattle Jar



Code	Description
WCS/397	Stainless steel, 203mm x 50mm

Measuring Jug




Code	Description
WCS/395	1 litre capacity

Petri Dish



Code	Description
WCS/002149	100mm, stainless steel, with lid

Gallipot



Code	Description
WCS/396	Stainless steel, 63mm



Wellcare surgicals . wellcaresurgicals@live.co.uk T. +44 7405903710 wellcaresurgicals.co.uk



Based in
the City of
Southend